

# SMARCAD1 Rabbit pAb

货号: **AYP15212**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:100
理论分子量	68kDa/117kDa
实测分子量	68kDa/117kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Chromosome,Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SMARCAD1 (NP_001121901.1).
序列	MNLFNLDLDRFRFEKRNKIEEAPQPSQPGPSSPISLSAEEENAEGEVSRANTPDSITEKTEDSSVPETPDNERKASISYFKNQRGIQYIDLSSDSEDVVSPNCSTVQEKTFNKDVTIIVSEPEDEESQGLPTMARRNDDISELEDLSELEDLKDACLQTLKELFPQRSNDLLKLIESTSTMDGAIAAALLMFGDAGGGPRKRKLSSESPEYEEDEFNDDQSIKKTRLDHGEESNESAESSNWEKQESIVLKLQKEFPNFDKQELREVLKEHEWMYTEALESKVFQEDQDMQY

## 靶点信息

研究背景	This gene encodes a member of the SNF subfamily of helicase proteins. The encoded protein plays a critical role in the restoration of heterochromatin organization and propagation of epigenetic patterns following DNA replication by mediating histone H3/H4 deacetylation. Mutations in this gene are associated with dermatoglyphia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
------	---

基因ID	56916
基因名	SMARCAD1
Swiss	Q9H4L7
别名	SMARCAD1;ADERM;BASNS;ETL1;HEL1

## 产品验证



Western blot analysis of SMARCAD1 expressed in human esophagus cancer, mouse kidney, rat liver using SMARCAD1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)