

IL11RA Rabbit pAb

货号: **AYP15241**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	41kDa/45kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse heart,Mouse skeletal muscle,Rat liver
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 30-290 of human IL11R A (NP_001136256.1).
序列	WGPPGVQYGQGRSVKLCPCGVTAGDPVSWFRDGEKLLQGPDSGLGHELVLQAADSTDEGTICQTLDGALGGTVTL QLGYPPARPVVSCQAADYENFCTWSPSQISGLPTRYLTSYRKKTVLGADSQRRSPSTGPWPCPDPLGAARCVVHGAE FWSQYRINVTEVNPLGASTRLLDVSLQSIIRPDPPQGLRVESVPGYPRRLRASWTYPASWPCQPHLLKFRLLQYRPAQHPA WSTVEPAGLEEVITDAVAGLPHA

靶点信息

研究背景	Interleukin 11 is a stromal cell-derived cytokine that belongs to a family of pleiotropic and redundant cytokines that use the gp130 transducing subunit in their high affinity receptors. This gene encodes the IL-11 receptor, which is a member of the hematopoietic cytokine receptor family. This particular receptor is very similar to ciliary neurotrophic factor, since both contain an extracellular region with a 2-domain structure composed of an immunoglobulin-like domain and a cytokine receptor-like domain. Multiple alternatively spliced transcript variants have been found for this gene.
------	--

基因ID	3590
基因名	IL11RA
Swiss	Q14626
别名	IL11RA;CRSDA

产品验证



Western blot analysis of IL11RA expressed in Mouse kidney,Rat liver,Mouse heart using IL11RA Rabbit p Ab at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s

实验步骤

访问官网浏览详情: www.ablybio.cn