

PCGF3 Rabbit pAb

货号: **AYP15291**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	28kDa
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse heart
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-242 of human PCGF3 (NP_006306.2).
序列	MLTRKIKLWDINAHITCRLCSGYLIDATTVTECLHTFCRSCLVKYLEENNTCPTCRIVIHQSHPLQYIGHDRMQDIVYKLVPG LQEAEMRKQREFYHKLGMVPGDIKGETCSAKQLDSDHRNGETKADDSSNKEAAEEKPEEDNDYHRSDEQVSICLECNS SKLRGLKRKWIRCSAQATVLHLKFKIAKLNLSFFNELDILCNEEILGKDHTLKFVVVTRWRFKKAPLLLHYRPMDDL

靶点信息

研究背景	The protein encoded by this gene contains a C3HC4 type RING finger, which is a motif known to be involved in protein-protein interactions. The specific function of this protein has not yet been determined.
基因ID	10336

基因名	PCGF3
Swiss	Q3KNV8
别名	PCGF3;DONG1;RNF3;RNF3A

产品验证



Western blot analysis of PCGF3 expressed in Mouse kidney, Mouse heart using PCGF3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn