

# STRN Rabbit pAb

货号: **AYP15298**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200
理论分子量	80kDa/86kDa
实测分子量	120kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,NIH/3T3,Mouse lung,Rat brain
细胞定位	Cell projection,Cytoplasm,Membrane,Peripheral membrane protein,dendritic spine
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 581-780 of human STRN (NP_003153.2).
序列	HQRLSCSADGTLRLWNTTEVAPALSVFNDDTKELGIPASVDLVSSDPSHMOVASFSGYTSIFNMETQQRILTLESNVDTTA NSSCQINRVISHPTLPISITAHEDRHIKFYDNNKGKLIHSMVAHLEAVTSLAVDPNGLYLMSGSHDCSIRLWNLESKTCIQEF TAHRKKFEESIHDFVAFHPSKCYIASAGADALAKVVF

## 靶点信息

研究背景	Calmodulin-binding protein which may function as scaffolding or signaling protein and may play a role in dendritic Ca <sup>2+</sup> signaling.
基因ID	6801

基因名	STRN
Swiss	O43815
别名	STRN;PPP2R6A;SG2NA;STRN1;striatin

## 产品验证



Western blot analysis of STRN expressed in Mouse lung,Rat brain,A-549 using STRN Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)