

— ABLYBIO, Help Your Research



STRN Rabbit pAb

货号: **AYP15298**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB IHC |
| 推荐浓度 | WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 |
| 理论分子量 | 80kDa/86kDa |
| 实测分子量 | 120kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | A-549,NIH/3T3,Mouse lung,Rat brain |
| 细胞定位 | Cell projection,Cytoplasm,Membrane,Peripheral membrane protein,dendritic spine |
| 纯化 | Affinity purification |


抗原信息

| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 581-780 of human STRN (NP_003153.2). |
|------|--|

靶点信息

| | |
|-------|--|
| 研究背景 | Calmodulin-binding protein which may function as scaffolding or signaling protein and may play a role in dendritic Ca ²⁺ signaling. |
| 基因ID | 6801 |
| 基因名 | STRN |
| Swiss | O43815 (https://www.uniprot.org/uniprotkb/O43815/entry) |
| 别名 | STRN,PPP2R6A,SG2NA,STRN1,striatin,STRN Rabbit pAb |

产品验证



Western blot analysis of STRN expressed in Mouse lung,Rat brain,A-549 using STRN Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)