

# PHF11 Rabbit pAb

货号: **AYP15301**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:100 - 1:500
理论分子量	33kDa/37kDa
实测分子量	38KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse liver,Mouse spleen
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of human PHF11 (NP_001035533.1).
序列	MAQASPPRPERVLGASSPEARPAQEALLLPTGVFQVAEKMEKRTCALCPKDVEYNVLYFAQSENIAAHENCLLYSSGLVECE DQDPLNPDRSFDVESVKKEIQRGRKLCCKFCHKRATVGCCLKNCNKNYHFFCAKKDDAVPQSDGVRGIYKLLCQQA QFPIIAQSAKFSGVKRKRGRKKPLSGNHVQPPEMTCNTFIRQVKEEHGRHTDATVKVPFLKKCKEAGLLNYLLEEILDKVV SIPKLMDETTSESDYEEIGSALFDCRLFEDTFVNFQAAIEKKIHASQQR

## 靶点信息

研究背景	This gene encodes a protein containing a PHD (plant homeodomain) type zinc finger. This gene has been identified in some studies as a candidate gene for asthma. Naturally-occurring readthrough transcription may occur from the upstream SETDB2 (SET domain bifurcated 2) gene to this locus. Alternative splicing results in multiple transcript variants.
------	---

基因ID	51131
基因名	PHF11
Swiss	Q9UIL8
别名	PHF11;APY;BCAP;IGEL;IGER;IGHER;NY-REN-34;NYREN34

## 产品验证



Western blot analysis of PHF11 expressed in Jurkat, Mouse liver, Mouse spleen using PHF11 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)