

— ABLYBIO, Help Your Research



PHF11 Rabbit pAb

货号: **AYP15301**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	33kDa/37kDa
实测分子量	38KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse liver,Mouse spleen
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of human PHF11 (NP_001035533.1).
------	--

靶点信息

研究背景	This gene encodes a protein containing a PHD (plant homeodomain) type zinc finger. This gene has been identified in some studies as a candidate gene for asthma. Naturally-occurring readthrough transcription may occur from the upstream SETDB2 (SET domain bifurcated 2) gene to this locus. Alternative splicing results in multiple transcript variants.
基因ID	51131
基因名	PHF11
Swiss	Q9UIL8 (https://www.uniprot.org/uniprotkb/Q9UIL8/entry)
别名	PHF11,APY,BCAP,IGEL,IGER,IGHER,NY-REN-34,NYREN34,PHF11 Rabbit pAb,BRCA1 C-terminus-associated protein,Renal carcinoma antigen NY-REN-34

产品验证



Western blot analysis of PHF11 expressed in Jurkat, Mouse liver, Mouse spleen using PHF11 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)