

— ABLYBIO, Help Your Research



# FTSJ1 Rabbit pAb

货号: **AYP15308**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	35kDa/36kDa
实测分子量	36KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse ovary,Rat thymus,Rat uterus
细胞定位	Cytoplasm
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 170-329 of human FTSJ 1 (NP_036412.1).
------	--

## 靶点信息

研究背景	This gene encodes a member of the methyltransferase superfamily. The encoded protein localizes to the nucleolus, binds to S-adenosylmethionine, and may be involved in the processing and modification of ribosomal RNA. Mutations in this gene are associated with mental retardation. Alternative splicing results in multiple transcript variants.
基因ID	24140
基因名	FTSJ1
Swiss	Q9UET6 ( <a href="https://www.uniprot.org/uniprotkb/Q9UET6/entry">https://www.uniprot.org/uniprotkb/Q9UET6/entry</a> )
别名	FTSJ1,CDLIV,JM23,MRX44,MRX9,SPB1,TRMT7,FTSJ1 Rabbit pAb,2'-O-ribose RNA methyltransferase TRM7 homolog,Protein ftsj homolog 1

## 产品验证



Western blot analysis of FTSJ1 expressed in HeLa, Mouse ovary, Rat thymus using FTSJ1 Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))