

PEX11B Rabbit pAb

货号: **AYP15369**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	28kDa
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse kidney,Mouse brain,Rat testis,Rat lung,Rat brain
细胞定位	cytosol,mitochondrion,nucleoplasm,peroxisome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-180 of human PEX11B (NP_003837.1).
序列	DNVLWAGKSGGLAPRVDQEKWAQRSFRYYLFLSLIMNLSRDAYEIRLLMEQESSACSRRLKSGGGVPGGSETGGLGGPG TPG

靶点信息

研究背景	The protein encoded by this gene facilitates peroxisomal proliferation and interacts with PEX19. The encoded protein is found in the peroxisomal membrane. Several transcript variants, some protein-coding and some not protein-coding, have been found for this gene.
基因ID	8799

基因名	PEX11B
Swiss	O96011
别名	PEX14B;PEX11-BETA;PEX11B

产品验证



Western blot analysis of PEX11B expressed in Mouse kidney,Rat testis,HepG2 using PEX11B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn