

— ABLYBIO, Help Your Research



# PEX11B Rabbit pAb

货号: **AYP15369**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	28kDa
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse kidney,Mouse brain,Rat testis,Rat lung,Rat brain
细胞定位	cytosol,mitochondrion,nucleoplasm,peroxisome
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-180 of human PEX11B (NP_003837.1).
------	--

## 靶点信息

研究背景	The protein encoded by this gene facilitates peroxisomal proliferation and interacts with PEX19. The encoded protein is found in the peroxisomal membrane. Several transcript variants, some protein-coding and some not protein-coding, have been found for this gene.
基因ID	8799
基因名	PEX11B
Swiss	O96011 ( <a href="https://www.uniprot.org/uniprotkb/O96011/entry">https://www.uniprot.org/uniprotkb/O96011/entry</a> )
别名	PEX14B,PEX11-BETA,PEX11B,PEX11B Rabbit pAb,Peroxin-11B,Peroxisomal biogenesis factor 11B,Protein PEX11 homolog beta

## 产品验证



Western blot analysis of PEX11B expressed in Mouse kidney,Rat testis,HepG2 using PEX11B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))