

CDY1B Rabbit pAb

货号: **AYP15394**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	60kDa/62kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse testis,Rat testis
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 60-320 of human CDY1 B (NP_001003895.1).
序列	EKQKCLTWTTSRIFSNARRRSTSRSTKANYSKNSPKTPVTDKHRSKNRKLFASKNVRRAASILSDTKNMEIINSTIETL APDSPFDHKT VSGFQKLEKLDPIAADQQDTVVFVKVTEGKLLRDPLSRPGAEQTGIQNKTIHPLMSQMSGSVTASMATG SATRKGIVVLIDPLAANGTTDMHTSVPRVKGGQRNITDDSRDQPFIKMHFTIRLTESASTYRDIVVKKEDGFTQIVLSTRS TEKNALNTEVIKEIVNA

靶点信息

研究背景	This gene encodes a protein containing a chromodomain and a histone acetyltransferase catalytic domain. Chromodomain proteins are components of heterochromatin-like complexes and can act as gene repressors. This protein is localized to the nucleus of late spermatids where histone hyperacetylation takes place. Histone hyperacetylation is thought to facilitate the transition in which protamines replace histones as the major DNA-packaging protein. The human chromosome Y has two identical copies of this gene within a palindromic region; this record represents the more centromeric copy. Chromosome Y also contains a pair of closely related genes in another more telomeric palindrome as well as several related pseudogenes. Two protein isoforms are encoded by transcript variants of this gene. Additional transcript variants have been described, but their full-length nature has not been determined.
基因ID	253175
基因名	
Swiss	Q9Y6F8
别名	CDY1B;CDY

产品验证



Western blot analysis of CDY1B expressed in Mouse testis, Rat testis using CDY1B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn