

— ABLYBIO, Help Your Research



# SUGP2 Rabbit pAb

货号: **AYP15417**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:200 - 1:2000
理论分子量	115kDa/116kDa/120kDa
实测分子量	131kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,293T,LO2,Jurkat,HeLa,Mouse thymus
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human SUG P2 (NP_055699.2).
------	--

## 靶点信息

研究背景	This gene encodes a member of the arginine/serine-rich family of splicing factors. The encoded protein functions in mRNA processing. Alternatively spliced transcript variants have been described.
------	---

基因ID	10147
基因名	SUGP2
Swiss	Q8IX01 ( <a href="https://www.uniprot.org/uniprotkb/Q8IX01/entry">https://www.uniprot.org/uniprotkb/Q8IX01/entry</a> )
别名	SUGP2,SFRS14,SUGP2 Rabbit pAb,Arginine/serine-rich-splicing factor 14,Splicing factor,arginine/serine-rich 14,KIAA0365

## 产品验证



Western blot analysis of SUGP2 expressed in Mouse thymus,U-251MG,293T using SUGP2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)