

— ABLYBIO, Help Your Research



PTP4A2 Rabbit pAb

货号: **AYP15422**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:200 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	9kDa/15kDa/16kDa/19kDa
实测分子量	19kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,Mouse skeletal muscle,Mouse brain,Mouse kidney,Mouse liver,Rat brain
细胞定位	Cell membrane,Cytoplasm,Early endosome
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 30-100 of human PTP4A 2 (NP_536316.1).
------	--

靶点信息

研究背景	The protein encoded by this gene belongs to a small class of the protein tyrosine phosphatase (PTP) family. PTPs are cell signaling molecules that play regulatory roles in a variety of cellular processes. PTPs in this class contain a protein tyrosine phosphatase catalytic domain and a characteristic C-terminal prenylation motif. This PTP has been shown to primarily associate with plasmic and endosomal membrane through its C-terminal prenylation. This PTP was found to interact with the beta-subunit of Rab geranylgeranyltransferase II (beta GGT II), and thus may function as a regulator of GGT II activity. Overexpression of this gene in mammalian cells conferred a transformed phenotype, which suggested its role in tumorigenesis. Alternatively spliced transcript variants have been described. Related pseudogenes exist on chromosomes 11, 12 and 17.
基因ID	8073
基因名	PTP4A2
Swiss	Q12974 (https://www.uniprot.org/uniprotkb/Q12974/entry)
别名	PTP4A2,HH13,HH7-2,HU-PP-1,OV-1,PRL-2,PRL2,PTP4A,PTPCAAX2,ptp-IV1a,ptp-IV1b,PTP4A2 Rabbit pAb,PTP(CAAXII),Protein-tyrosine phosphatase 4a2,Protein-tyrosine phosphatase of regenerating liver 2

产品验证



Western blot analysis of PTP4A2 expressed in Mouse skeletal muscle,Rat brain,MCF7 using PTP4A2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)