

# ZNF677 Rabbit pAb

货号: **AYP15444**

## 产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	67kDa
实测分子量	68kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Rat testis
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 50-200 of human ZNF677 (NP_872415.1).
序列	DNIPPEDDISVGFTSKGLSPKENNKEELYHLVILERKESHGINNFDLKEVWENMPKFDLSLWDYDVKNYKGMPLTCNKNLTH RKDQQHNKSSIHFSLKQSVSIRDSAHQYFIHDKPFIRNLLKLNIRYAGNKYVKCFENKIGLSLQAQLA

## 靶点信息

研究背景	May be involved in transcriptional regulation.
基因ID	342926
基因名	ZNF677

Swiss	Q86XU0
别名	ZNF677

## 产品验证



Western blot analysis of ZNF677 expressed in Rat testis using ZNF677 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)