

— ABLYBIO, Help Your Research



IL31 Rabbit pAb

货号: **AYP15445**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	18kDa
实测分子量	18kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HT-29,Mouse spleen,Mouse kidney,Mouse thymus,Rat testis
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 24-164 of human IL31 (NP_001014358.1).
------	--

靶点信息

研究背景	
------	--

基因ID	386653
基因名	IL31
Swiss	Q6EBC2 (https://www.uniprot.org/uniprotkb/Q6EBC2/entry)
别名	IL31,IL-31,IL31 Rabbit pAb

产品验证



Western blot analysis of IL31 expressed in Mouse spleen,Rat testis,293T using IL31 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)