

# ACTL8 Rabbit pAb

货号: **AYP15450**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	41kDa
实测分子量	41-50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2,Mouse testis,Mouse pancreas,Mouse brain
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-366 of human ACTL8 (NP_110439.2).
序列	MAARTVIIDHGSGFLKAGTAGWNEPQMVFVNIVNYLPCKENPGPSYARRRVSLGIDICHPDTFSYPIERGRILNWEGVQYL WSFVLENHRREQEVPPVITETPLREPADRKKMLEILFELLHVPVSVLLADQLQMSLYASGLLTGVVVDVSGYGLTRVQPFHQG RPLPASGKTLEFAGQDLSAYLLKSLFKEDCDRRCLFQLETVAVTQMNCYVPQNLGEALDFRERQQSALDESNTYQLPDG SRVELTPMQRVAPEMFFSPQVFEQPGPSIPRAIVESVESCEISLRPLLVSHVMACGGNTLYPGFTKRLFRELMGDHVSSTKA TVWEGSNRNFVWLGASVVAHLSTYQSEWMSREEYGEHMRM

## 靶点信息

研究背景	Involved in epithelial cell differentiation. Predicted to be located in cytoplasm. Predicted to be part of dyn actin complex.
基因ID	81569

基因名	ACTL8
Swiss	Q9H568
别名	ACTL8;CT57

## 产品验证



Western blot analysis of ACTL8 expressed in Mouse testis, HeLa, HepG2 using ACTL8 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)