

# LYPD3 Rabbit pAb

货号: **AYP15455**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	35kDa
实测分子量	36kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3, MCF7, HT-29
细胞定位	Cell membrane, GPI-anchor, Lipid-anchor
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 31-326 of human LYPD3 (NP_055215.2).
序列	LECYSCVQKADDGCSPNKMKTVKCAPGVDVCTEAVGAVETIHGQFSLAVRGCGLPGKNDRGLDLHGLLAFIQLQOCA QDRCNAKLNLTSRALDPAGNESAYPPNGVECYSCVGLSREACQGTSPVWSCYNASDHVYKGCDFGNVTLTAANVTVSL PVRGCVQDEFCTRDGVTGPGFTLSGCCQGSRCNSDLRNKTYFSPRIPLVRLPPPEPTTVASTTSVTTSTSTAPVRPTSTT KMPAPTSQTPRQGVHEASRDEEPRLTGGAAGHQDRSNGQYPAKGGPQQPHNKGC

## 靶点信息

研究背景	Supports cell migration. May be involved in urothelial cell-matrix interactions. May be involved in tumor progression.
基因ID	27076

基因名	LYPD3
Swiss	O95274
别名	LYPD3;C4.4A

## 产品验证



Western blot analysis of LYPD3 expressed in BxPC-3, MCF7, HT-29 using LYPD3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)