

— ABLYBIO, Help Your Research



KIRREL Rabbit pAb

货号: **AYP15498**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	72kDa/83kDa/85kDa
实测分子量	83kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse lung
细胞定位	Cell membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 17-200 of human KIRREL (NP_060710.3).
------	---

靶点信息

研究背景	
------	--

基因ID	55243
基因名	KIRREL
Swiss	Q96J84 (https://www.uniprot.org/uniprotkb/Q96J84/entry)
别名	KIRREL, NEPH1, KIRREL Rabbit pAb, Kin of irregular chiasm-like protein 1, Nephrin-like protein 1

产品验证



Western blot analysis of KIRREL expressed in Mouse kidney, Mouse lung using KIRREL Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)