


基因名	FXR2
Swiss	P51116
别名	FXR2;FMR1L2;FXR2P

产品验证



Western blot analysis of FXR2 expressed in Mouse testis,Rat liver,MCF7 using FXR2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn