

— ABLYBIO, Help Your Research



PTN Rabbit pAb

货号: **AYP15509**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	18kDa
实测分子量	19kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,MCF7,NIH/3T3,Mouse liver
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human PTN (NP_002816.1).
------	--

靶点信息

研究背景	
------	--

基因ID	5764
基因名	PTN
Swiss	P21246 (https://www.uniprot.org/uniprotkb/P21246/entry)
别名	PTN,HARP,HB-GAM,HBBM,HBGF-8,HBGF8,HBNF,HBNF-1,NEGF1,OSF-1,PTN Rabbit pAb,Heparin-binding bra in mitogen,Heparin-binding growth factor 8,Heparin-binding growth-associated molecule,Heparin-binding neurite outgrowth-promoting factor,Heparin-binding neurite outgrowth-promoting factor 1

产品验证



Western blot analysis of PTN expressed in Mouse liver,U-87MG,MCF7 using PTN Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)