

— ABLYBIO, Help Your Research



# ARHGAP25 Rabbit pAb

货号: **AYP15524**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	51kDa/69kDa/72kDa/73kDa
实测分子量	73kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,A-549,Mouse lung,Rat spleen
细胞定位	phagocytic cup
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-370 of human ARHGA P25 (NP_001159748.1).
------	--

## 靶点信息

研究背景	
------	--

基因ID	9938
基因名	ARHGAP25
Swiss	P42331 ( <a href="https://www.uniprot.org/uniprotkb/P42331/entry">https://www.uniprot.org/uniprotkb/P42331/entry</a> )
别名	ARHGAP25,HEL-S-308,KAIA0053,ARHGAP25 Rabbit pAb,Rho-type GTPase-activating protein 25,KIAA0053

## 产品验证



Western blot analysis of ARHGAP25 expressed in Mouse lung,Rat spleen,Jurkat using ARHGAP25 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)