

# ACTRT3 Rabbit pAb

货号: **AYP15545**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	41kDa
实测分子量	41kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HT-29,293T,Mouse testis,Mouse brain,Mouse kidney,Rat testis,Rat brain,Rat kidney
细胞定位	Cytoplasm,Nucleus,cytoskeleton
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human ACTRT3 (NP_115876.3).
序列	MNHCQLPVVIDNGSGMIKAGVAGCREPQFIYPNIIGRAKQSRRAQGGLELCVGDQAQDWRSSLFISYPVERGLITSWED MEIMWKHIYDYNLKLKPCDGPVLITEPALNPLANRQQITEMFFEHLGVPA

## 靶点信息

研究背景	Predicted to be located in male germ cell nucleus.
基因ID	84517
基因名	ACTRT3

<b>Swiss</b>	Q9BYD9
别名	ACTRT3;ARP-T3;ARPM1

## 产品验证



Western blot analysis of ACTRT3 expressed in Mouse testis,Rat testis,U-87MG using ACTRT3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)