

SAAL1 Rabbit pAb

货号: **AYP15546**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	53kDa
实测分子量	54-58kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,LO2,HeLa,Jurkat,Mouse lung,Mouse pancreas,Rat kidney
细胞定位	nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 175-474 of human SAAL 1 (NP_612430.2).
序列	ERIQEHPAIYDSICFIMSSSTNVDLLVKVGEVVDKLFDLDEKLMLEWVRNGAAQPLDQPQESEEQPVFRLVPCILEAAKQV RSENPEWLDVYMHLQLLTTVDDGIQAIVHCPDTGKDIWNLLFDLVCHEFCQSDDPPIILQEQTVLASVFSVLSAIYASQT EQEYLKIEKVDLPLIDSLIRVLQNMEQCQKPKPNSAESNTEETKRDTLTQDDFHLKILKDILCEFLSNIFQALTKETVAQGVKE GQLSKQKCSSAFQNLPPFYSPVVEDFIKILREVDKALADDLEKNFPSLKVQT

靶点信息

研究背景	Plays a role in promoting the proliferation of synovial fibroblasts in response to proinflammatory stimuli.
基因ID	113174
基因名	SAAL1

Swiss	Q96ER3
别名	SAAL1;SPACIA1

产品验证



Western blot analysis of SAAL1 expressed in Mouse lung,Rat kidney,293T using SAAL1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn