

UBAC1 Rabbit pAb

货号: **AYP15558**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	45kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse liver,Mouse lung,Rat testis
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 176-405 of human UBA C1 (NP_057256.2).
序列	LFKKANAMLDEDEDERVDAAALRQLTEMGFPENRATKALQLNHMSVPQAMEWLIHAEDPTIDTPLPGQAPPEAEGATAA ASEAAAAGASATDEEARDELTEIFKKIRRKREFRADARAVISLMEMGFDEKEVIDALRVNNNQNAACEWLLGDRKPSPEEL DKGIDPDSPLFQAILDNPVVQLGLTNPKTLAFEDMLENPLNSTQWMNDPETGPVMLQISRIFQTLNRT

靶点信息

研究背景	Non-catalytic subunit of the KPC complex that acts as E3 ubiquitin-protein ligase. Required for poly-ubiquitination and proteasome-mediated degradation of CDKN1B during G1 phase of the cell cycle.
基因ID	10422

基因名	UBAC1
Swiss	Q9BSL1
别名	UBAC1;GBDR1;KPC2;UBADC1

产品验证



Western blot analysis of UBAC1 expressed in Mouse liver,Rat testis,Jurkat using UBAC1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn