

— ABLYBIO, Help Your Research



# SNRNP40 Rabbit pAb

货号: **AYP15669**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	
实测分子量	40KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa,Mouse brain,Mouse testis
细胞定位	cytosol,nuclear speck,nucleoplasm,nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-357 of human SNRNP 40 (NP_004805.2).
------	--

## 靶点信息

研究背景	This gene encodes a component of the U5 small nuclear ribonucleoprotein (snRNP) particle. The U5 snRNP is part of the spliceosome, a multiprotein complex that catalyzes the removal of introns from pre-messenger RNAs. [provided by RefSeq, Jul 2008]
基因ID	9410
基因名	SNRNP40
Swiss	Q96DI7 ( <a href="https://www.uniprot.org/uniprotkb/Q96DI7/entry">https://www.uniprot.org/uniprotkb/Q96DI7/entry</a> )
别名	40K,HPRP8BP,PRP8BP,PRPF8BP,SPF38,WDR57,SNRNP40 Rabbit pAb,SNRNP40,38 kDa-splicing factor,Prp8-binding protein,U5 snRNP-specific 40 kDa protein,WD repeat-containing protein 57,SFP38

## 产品验证



Western blot analysis of SNRNP40 expressed in Mouse brain,293T,HeLa using SNRNP40 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)