


基因ID	9410
基因名	SNRNP40
Swiss	Q96DI7
别名	40K; HPRP8BP; PRP8BP; PRPF8BP; SPF38; WDR57

产品验证



Western blot analysis of SNRNP40 expressed in Mouse brain, 293T, HeLa using SNRNP40 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn