

— ABLYBIO, Help Your Research



# PDHA2 Rabbit pAb

货号: **AYP15703**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	42kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NCI-H460,SKOV3,A-549,Mouse testis,Mouse ovary
细胞定位	Mitochondrion matrix
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 119-388 of human PDH A2 (NP_005381.1).
------	--

## 靶点信息

研究背景	The pyruvate dehydrogenase complex catalyzes the overall conversion of pyruvate to acetyl-CoA and CO <sub>2</sub> , and thereby links the glycolytic pathway to the tricarboxylic cycle.
基因ID	5161
基因名	PDHA2
Swiss	P29803 ( <a href="https://www.uniprot.org/uniprotkb/P29803/entry">https://www.uniprot.org/uniprotkb/P29803/entry</a> )
别名	PDHA2,PDHAL,PDHA2 Rabbit pAb,PDHE1-A type II

## 产品验证



Western blot analysis of PDHA2 expressed in Mouse testis, NCI-H460, SKOV3 using PDHA2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))