

# ZFPM2 Rabbit pAb

货号: **AYP15728**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	113kDa/128kDa
实测分子量	150kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-1080,K-562,Mouse liver
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 982-1151 of human ZFP M2 (NP_036214.2).
序列	SPYYGIKPSDYISGSLVIHNTDIEQSRNAENESPKGQASSNGCAALKKDSLPLLPKNRGMVIVNGGLKQDERPAANPQQEN ISQNPQHEDDHKSPSWISENPLAANENVSPGIPSAEEQLSSIAKGVNGSSQAPTSQKRCRLCDIQFNLSNFITHKKFYCSS HAAEHVK

## 靶点信息

研究背景	The zinc finger protein encoded by this gene is a widely expressed member of the FOG family of transcription factors. The family members modulate the activity of GATA family proteins, which are important regulators of hematopoiesis and cardiogenesis in mammals. It has been demonstrated that the protein can both activate and down-regulate expression of GATA-target genes, suggesting different modulation in different promoter contexts. A related mRNA suggests an alternatively spliced product but this information is not yet fully supported by the sequence.
------	--

基因ID	23414
基因名	ZFPM2
Swiss	Q8WW38
别名	ZFPM2;DIH3;FOG2;SRXY9;ZC2HC11B;ZNF89B;hFOG-2

## 产品验证



Western blot analysis of ZFPM2 expressed in HT-1080, K-562, Mouse liver using ZFPM2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)