

— ABLYBIO, Help Your Research



# FBN2 Rabbit pAb

货号: **AYP15803**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	157kDa/314kDa
实测分子量	290kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,Mouse brain
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 2733-2912 of human FB N2 (NP_001990.2).
------	---

## 靶点信息

研究背景	The protein encoded by this gene is a component of connective tissue microfibrils and may be involved in elastic fiber assembly. Mutations in this gene cause congenital contractural arachnodactyly.
------	---

基因ID	2201
基因名	FBN2
Swiss	P35556 ( <a href="https://www.uniprot.org/uniprotkb/P35556/entry">https://www.uniprot.org/uniprotkb/P35556/entry</a> )
别名	FBN2,CCA,DA9,EOMD,FBN2 Rabbit pAb

## 产品验证



Western blot analysis of FBN2 expressed in LO2, Mouse brain using FBN2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)