

FBN2 Rabbit pAb

货号: **AYP15803**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:2000 |
| 理论分子量 | 157kDa/314kDa |
| 实测分子量 | 290kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | LO2,Mouse brain |
| 细胞定位 | Secreted,extracellular matrix,extracellular space |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 2733-2912 of human FB N2 (NP_001990.2). |
| 序列 | CVSGMGFNKGQYLSLDTEVDEENALSPEACYECKINGYSKKDSRQKRSIHEPDPTAVEQISLESVDMDSPVNMKFNLSHLGSKEHILELRPAIQPLNNHIRYVISQGNDDSVFRIHQRNGLSYLHTAKKKLMPGTYTLEITSIPLYKKKELKKLEESNEDDYLLGELGEALRMRLQIQLY |

靶点信息

| | |
|------|---|
| 研究背景 | The protein encoded by this gene is a component of connective tissue microfibrils and may be involved in elastic fiber assembly. Mutations in this gene cause congenital contractural arachnodactyly. |
| 基因ID | 2201 |

| | |
|-------|-------------------|
| 基因名 | FBN2 |
| Swiss | P35556 |
| 别名 | FBN2;CCA;DA9;EOMD |

产品验证



Western blot analysis of FBN2 expressed in LO2, Mouse brain using FBN2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn