

— ABLYBIO, Help Your Research



# SNRPB2 Rabbit pAb

货号: **AYP15871**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	25kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,U-87MG,HT-1080,LO2,293T,Mouse liver
细胞定位	Nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human SNRPB 2 (NP_003083.1).
------	---

## 靶点信息

研究背景	The protein encoded by this gene associates with stem loop IV of U2 small nuclear ribonucleoprotein (U2 snRNP) in the presence of snRNP-A'. The encoded protein may play a role in pre-mRNA splicing. Autoantibodies from patients with systemic lupus erythematosus frequently recognize epitopes on the encoded protein. Two transcript variants encoding the same protein have been found for this gene.
基因ID	6629
基因名	SNRPB2
Swiss	P08579 ( <a href="https://www.uniprot.org/uniprotkb/P08579/entry">https://www.uniprot.org/uniprotkb/P08579/entry</a> )
别名	SNRPB2,Msl1,U2B,SNRPB2 Rabbit pAb

## 产品验证



Western blot analysis of SNRPB2 expressed in Mouse liver, HeLa, U-87MG using SNRPB2 Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)