

# SRFBP1 Rabbit pAb

货号: **AYP15878**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	48kDa
实测分子量	48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart, Mouse skeletal muscle, Mouse liver, Mouse kidney
细胞定位	Cytoplasm, perinuclear region
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SRFBP1 (NP_689759.2).
序列	MAQPGTLNLNNEVVKMRKEVKRIRVLVIRKLVRSVGRLLKSKKGTEDALLKNQRRARQLLEEIHAMKELKPDIVTKSALGDDI NFEKIFKKPDSTATERAIARLAVHPLLKKKIDVLKAAVQAFKEARQNVAEVSSKNASEDNHSENTLYSNDNGSNLQREATV ISEQVKETKILAKKPIHNSKEKIAKMEHGPKAVTIANSPSPKSEKDSVVSLESQKTPADPKLKTLSQTKKNKGSDDSSLSGNS DGGEEFCHEEKEYFDDSTEERFYKQSSMSSESDSGDDFFIGKVRRTKKESSC

## 靶点信息

研究背景	May be involved in regulating transcriptional activation of cardiac genes during the aging process. May play a role in biosynthesis and/or processing of SLC2A4 in adipose cells (By similarity).
基因ID	153443

基因名	SRFBP1
Swiss	Q8NEF9
别名	SRFBP1;BUD22;P49;Rib1;STRAP;p49/STRAP

## 产品验证



Western blot analysis of SRFBP1 expressed in Mouse heart, Mouse skeletal muscle, Mouse liver using SRFBP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/prot eins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)