

# CENPJ Rabbit pAb

货号: **AYP15898**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	124kDa/153kDa
实测分子量	153kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji
细胞定位	Cytoplasm,centriole,centrosome,cytoskeleton,microtubule organizing center
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1069-1338 of human CE NPJ (NP_060921.3).
序列	DKLANTSVRFQNSQISSGTQVEKYKKNYLPMQGNPPRRSKSAPPRDLGNLDKGQAASPREPLEPLNFPDPEYKEEEDQD IQGEISHPDGKVEKVYKNGCRVILFPNGTRKEVSADGKTITVFFNGDVKQVMPDQRVIYAAAQTTHHTYPEGLEVLHFS SGQIEKHYPDGRKEITFPDQTVKNLFPDQGEESIFPDGTIVRVQRDGNKLIENNGQRELHTAQFKRREYPDGTVKTVYAN GHQETKYRSGRIRVKDEGNVLMDELT

## 靶点信息

研究背景	This gene encodes a protein that belongs to the centromere protein family. During cell division, this protein plays a structural role in the maintenance of centrosome integrity and normal spindle morphology, and it is involved in microtubule disassembly at the centrosome. This protein can function as a transcriptional coactivator in the Stat5 signaling pathway, and also as a coactivator of NF-kappaB-mediated transcription, likely via its interaction with the coactivator p300/CREB-binding protein. Mutations in this gene are associated with primary autosomal recessive microcephaly, a disorder characterized by severely reduced brain size and mental retardation. Alternatively spliced transcript variants have been found for this gene.
基因ID	55835
基因名	CENPJ
Swiss	Q9HC77
别名	CENPJ;BM032;CENP-J;CPAP;LAP;LIP1;MCPH6;SASS4;SCKL4;Sas-4

## 产品验证



Western blot analysis of CENPJ expressed in Raji using CENPJ Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)