

# MIB1/DIP1 Rabbit pAb

货号: **AYP15925**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:100 - 1:200
理论分子量	110kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A-549,HepG2,Mouse stomach,Mouse lung,Mouse kidney ,Rat testis
细胞定位	Cell membrane,Cytoplasm,centriolar satellite,centrosome,cytoskeleton,microtubule organizing center
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 560-860 of human MIB1 /DIP1 (NP_065825.1).
序列	DSEGDTPHDAISKRRDDILAVLLEAGADVTTITNNGFNALHHAALRGNPSSAMRVLLSKLPRPWIVDEKKDDGYTALHLAA LNNHVEVAELLVHQGNANLDIQNVNQQTALHLAVERQHTQIVRLLVRAGAKLDIQDKDGDTPLHEALRHHTLSQLRQLQD MQDVGKVDAAWEPKNTLIMGLGTQGAEKSAASIACFLAANGADLSIRNKKGQSPDLCPDNLCKALAKCHKKEKVS QVGSRSPTMISNDSETLEECMVCSDMKRDTLFGPCGHIATCSLCSRPVKKCLICKEQVQSR

## 靶点信息

研究背景	This gene encodes a protein containing multiple ankyrin repeats and RING finger domains that functions as an E3 ubiquitin ligase. The encoded protein positively regulates Notch signaling by ubiquitinating the Notch receptors, thereby facilitating their endocytosis. This protein may also promote the ubiquitination and degradation of death-associated protein kinase 1 (DAPK1).
------	--

基因ID	57534
基因名	MIB1
Swiss	Q86YT6
别名	MIB1;DIP-1;DIP1;LVNC7;MIB;ZZANK2;ZZZ6

## 产品验证



Western blot analysis of MIB1/DIP1 expressed in Mouse stomach,Rat testis,U-87MG using MIB1/DIP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)