

— ABLYBIO, Help Your Research



# HERC3 Rabbit pAb

货号: **AYP15934**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	39kDa/117kDa
实测分子量	117kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,U-87MG,A-431,Mouse lung,Mouse spinal cord
细胞定位	Cytoplasm,Cytoplasmic vesicle
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human HERC3 (NP_055421.1).
------	---

## 靶点信息

研究背景	This gene encodes a member the HERC ubiquitin ligase family. The encoded protein is located in the cytosol and binds ubiquitin via a HECT domain. Mutations in this gene have been associated with colorectal and gastric carcinomas. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	8916
基因名	HERC3
Swiss	Q15034 ( <a href="https://www.uniprot.org/uniprotkb/Q15034/entry">https://www.uniprot.org/uniprotkb/Q15034/entry</a> )
别名	HERC3,HERC3 Rabbit pAb,HECT domain and RCC1-like domain-containing protein 3,HECT-type E3 ubiquitin transferase HERC3,KIAA0032

## 产品验证



Western blot analysis of HERC3 expressed in Mouse lung,SH-SY5Y,U-87MG using HERC3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))