

— ABLYBIO, Help Your Research



UTY Rabbit pAb

货号: **AYP15935**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	118kDa/134kDa/137kDa/149kDa/154kDa
实测分子量	182kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse testis,Mouse liver,Rat kidney
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human UTY (NP_009056.3).
------	--

靶点信息

研究背景	This gene encodes a protein containing tetratricopeptide repeats which are thought to be involved in protein-protein interactions. The encoded protein is also a minor histocompatibility antigen which may induce graft rejection of male stem cell grafts. A large number of alternatively spliced transcripts have been observed for this gene, but the full length nature of some of these variants has not been determined.
基因ID	7404
基因名	UTY
Swiss	O14607 (https://www.uniprot.org/uniprotkb/O14607/entry)
别名	UTY,KDM6AL,KDM6C,UTY1,UTY Rabbit pAb,Ubiquitously-transcribed TPR protein on the Y chromosome,Ubiquitously-transcribed Y chromosome tetratricopeptide repeat protein,[histone H3]-trimethyl-L-lysine(27) demethylase UTY

产品验证



Western blot analysis of UTY expressed in Mouse kidney,Rat kidney,Mouse testis using UTY Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)