

ANAPC7 Rabbit pAb

货号: **AYP15947**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	
实测分子量	67KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	cytoplasm,cytosol,microtubule cytoskeleton,nucleoplasm,nucleus,spindle
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human ANAPC7 (NP_057322.2).
序列	HSNVRLLSLLLTMSNNPELFSPQKYQLLVYHADSLFHDKEYRNAVSKYTMALQKKALSKTSKVRPSTGNSASTPQS QCLPSEIEVKYKMAECYMLK

靶点信息

研究背景	This gene encodes a tetratricopeptide repeat containing component of the anaphase promoting complex/cyclosome (APC/C), a large E3 ubiquitin ligase that controls cell cycle progression by targeting a number of cell cycle regulators such as B-type cyclins for 26S proteasome-mediated degradation through ubiquitination. The encoded protein is required for proper protein ubiquitination function of APC/C and for the interaction of APC/C with certain transcription coactivators. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2008]
------	--

基因ID	51434
基因名	ANAPC7
Swiss	Q9UJX3
别名	APC7;ANAPC7

产品验证



Western blot analysis of ANAPC7 expressed in HepG2 using ANAPC7 Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn