

— ABLYBIO, Help Your Research



EMCN Rabbit pAb

货号: **AYP15967**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	18kDa/26kDa/27kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,NIH/3T3
细胞定位	Cell membrane,Membrane,Secreted,Single-pass type I membrane protein,Single-pass type I membrane p rotein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 19-135 of human EMCN (NP_057326.2).
------	---

靶点信息

研究背景	EMCN is a mucin-like sialoglycoprotein that interferes with the assembly of focal adhesion complexes and inhibits interaction between cells and the extracellular matrix (Kinoshita et al., 2001 [PubMed 11418125]).[supplied by OMIM, Mar 2008]
基因ID	51705
基因名	EMCN
Swiss	Q9ULC0 (https://www.uniprot.org/uniprotkb/Q9ULC0/entry)
别名	EMCN,EMCN2,MUC14,endomucin,EMCN Rabbit pAb,Endomucin-2,Gastric cancer antigen Ga34,Mucin-14

产品验证



Western blot analysis of EMCN expressed in MCF7,NIH/3T3 using EMCN Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)