

# RNF31 Rabbit pAb

货号: **AYP16061**

## 产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:100
理论分子量	57kDa/102kDa/119kDa
实测分子量	125kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,MCF7,BT-474,HeLa,PC-3
细胞定位	Cytoplasm
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human RNF31 (NP_060469.4).
序列	AAGACPEEIFSALQYSGTEVPLQWLRSELPYVLEMVAELAGQQDPGLGAFSCQEARRAWLDRHGNLDEAVEECVRTRRR KVQELQSLGFGPEEGSLQALFQHGGDVSRA TELQRQRLEPFRQRLWDSGPEPTPSWDGDPKQSLVRRLLAVYALPSW GRAELALSLLQETPRNYELGDVVEAVRHSQDRAFLRLLAQECA

## 靶点信息

研究背景	The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions. The encoded protein is the E3 ubiquitin-protein ligase component of the linear ubiquitin chain assembly complex. Two transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	55072
基因名	RNF31
Swiss	Q96EP0
别名	RNF31;HOIP;Paul;ZIBRA

## 产品验证



Western blot analysis of RNF31 expressed in Jurkat, MCF7, BT-474 using RNF31 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)