

# ANXA8L2 Rabbit pAb

货号: **AYP16259**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	30kDa/36kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HeLa,Mouse lung,Mouse ovary
细胞定位	
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-276 of human ANXA8L2 (NP_001265853.1).
序列	MAWWKAWIEQEGVTVKSSSHFNPDPAETLYKAMKGIGVGSQLLSHQAAAFAPSSALTSVSPWGQQGHLCCNPAGT NEQAIIDVLTKRSNTQRQQAQFKAQFGKDLTETLKSELGKFERLIVALMYPYRYEAKELHDAMKGSRDDVSSFVDPAL ALQDAQDLYAAGEKIRGTDEMKFITILCTRSATHLLRVKCTQNLHSYFAERLYYAMKGAGTRDGLRINIIVSRSEIDLNLKIC HFKKMYGKTLSSMIMEDTSGDYKNALLSLVGSDP

## 靶点信息

研究背景	This gene encodes a member of the annexin family of evolutionarily conserved Ca <sup>2+</sup> and phospholipid binding proteins. The encoded protein may function as an anticoagulant that indirectly inhibits the thromboplastin-specific complex. Overexpression of this gene has been associated with acute myelocytic leukemia. A highly similar duplicated copy of this gene is found in close proximity on the long arm of chromosome 10.
------	---

基因ID	728113
基因名	ANXA8L1
Swiss	Q5VT79
别名	ANXA8L1;ANXA8;ANXA8L2;VAC-beta;bA145E20.2

## 产品验证



Western blot analysis of ANXA8L2 expressed in Mouse lung,A-549,HeLa using ANXA8L2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)