

DVL2 Rabbit pAb

货号: **AYP16290**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:200
理论分子量	78kDa
实测分子量	78kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell membrane,Cytoplasm,Cytoplasmic side,Cytoplasmic vesicle,Peripheral membrane protein,cytosol
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-340 of human DVL2 (NP_004413.1).
序列	MAGSSTGGGGVGETKVIYHLDEEETPYLVKIPVPAERITLGDFKSVLQRPAGAKYFFKSMQDFGVVKEEISDDNARLPCF NGRVVSWLVSSDNPQPEMAPPVHEPRAELAPPAPLPLPPERTSGIGDSRPPSFHPNVSSHENLEPETETESVSLRRER PRRRDSSEHGAGGHRTGGPSRLERHLAGYESSTLMTSELESTSLGDSDEEDTMSRFSSTEQSSASRLLKRHHRRRKQ RPPRLERTSSFSSVTDSTMSLNIITVTLNMEKYNFLGISIVGQSNERGDGGIYIGSIMKGGAVAADGRIEPGDMLLQVNDM NFENMSNDDAVRVLRLDI

靶点信息

研究背景	This gene encodes a member of the dishevelled (dsh) protein family. The vertebrate dsh proteins have approximately 40% amino acid sequence similarity with Drosophila dsh. This gene encodes a 90-kD protein that undergoes posttranslational phosphorylation to form a 95-kD cytoplasmic protein, which may play a role in the signal transduction pathway mediated by multiple Wnt proteins. The mechanisms of dishevelled function in Wnt signaling are likely to be conserved among metazoans.
------	--

基因ID	1856
基因名	DVL2
Swiss	O14641
别名	DVL2

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn