

— ABLYBIO, Help Your Research



# DVL2 Rabbit pAb

货号: **AYP16290**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	<b>IF/ICC:</b> 1:50 - 1:200
理论分子量	78kDa
实测分子量	78kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell membrane,Cytoplasm,Cytoplasmic side,Cytoplasmic vesicle,Peripheral membrane protein,cytosol
纯化	Affinity purification

## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-340 of human DVL2 (NP_004413.1).
------	--

## 靶点信息

研究背景	This gene encodes a member of the dishevelled (dsh) protein family. The vertebrate dsh proteins have approximately 40% amino acid sequence similarity with Drosophila dsh. This gene encodes a 90-kD protein that undergoes posttranslational phosphorylation to form a 95-kD cytoplasmic protein, which may play a role in the signal transduction pathway mediated by multiple Wnt proteins. The mechanisms of dishevelled function in Wnt signaling are likely to be conserved among metazoans.
基因ID	1856
基因名	DVL2
Swiss	O14641 ( <a href="https://www.uniprot.org/uniprotkb/O14641/entry">https://www.uniprot.org/uniprotkb/O14641/entry</a> )
别名	DVL2,DVL2 Rabbit pAb,DSH homolog 2

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)