

— ABLYBIO, Help Your Research



# PPP3R2 Rabbit pAb

货号: **AYP16312**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	19kDa
实测分子量	16kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,H9C2,Mouse heart,Rat heart
细胞定位	cytosol,mitochondrion
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-173 of human PPP3R2 (NP_671709.1).
------	--

## 靶点信息

研究背景	Regulatory subunit of calcineurin, a calcium-dependent, calmodulin stimulated protein phosphatase. Confers calcium sensitivity.
------	---

基因ID	5535
基因名	PPP3R2
Swiss	Q96LZ3 ( <a href="https://www.uniprot.org/uniprotkb/Q96LZ3/entry">https://www.uniprot.org/uniprotkb/Q96LZ3/entry</a> )
别名	PPP3R2,PPP3RL,PPP3R2 Rabbit pAb,Calcineurin B-like protein,Calcineurin BII,PPP3R1-like,Protein phosphatase 2B regulatory subunit 2,Protein phosphatase 3 regulatory subunit B beta isoform,CBLP

## 产品验证



Western blot analysis of PPP3R2 expressed in Mouse heart,Rat heart,HeLa using PPP3R2 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blo cking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)