

CYP2C18 Rabbit pAb

货号: **AYP16344**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	48kDa/55kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	Endoplasmic reticulum membrane,Microsome membrane,Peripheral membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 231-490 of human CYP2C18 (NP_000763.1).
序列	NKIAENFAYIKSYVLERIKEHQESLDMNSARDFIDCFLIKMEQEKHNNQSEFTVESLIATVTDMFGAGTETTSTTLRYGLLLLL KYPEVTAKVQEEIECVVGRNRSRSPCMQDRSHMPYTDVAVVHEIQRYIDLLPTNLPHAVTCDVKFKNYLIPKGTIITSLTSVLHN DKEFPNPEMFDPGHFLDKSGNFKKSDFMPFSAGKRMCMGEGLARMELFLFTTILQNFNLKSQVDPKDIDITPIANAFGR VPPLYQLCFIPV

靶点信息

研究背景	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum but its specific substrate has not yet been determined. The gene is located within a cluster of cytochrome P450 genes on chromosome 10q24. An additional gene, CYP2C17, was once thought to exist; however, CYP2C17 is now considered an artefact based on a chimera of CYP2C18 and CYP2C19. Alternatively spliced transcripts encoding different isoforms have been found for this gene.
基因ID	1562
基因名	CYP2C18
Swiss	P33260
别名	CYP2C18;CPCI;CYP2C;CYP2C17;P450-6B/29C;P450IIC17

产品验证



Western blot analysis of CYP2C18 expressed in Mouse liver,Rat liver using CYP2C18 Rabbit pAb at 1:1000 . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn