

AEBP2 Rabbit pAb

货号: **AYP16354**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	33kDa/52kDa/54kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,293T,A-549,Mouse kidney,Mouse skeletal muscle,Rat heart
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 234-503 of human AEBP 2 (NP_694939.2).
序列	STISSGRSTPAMMNGQGSTTSSSKNIAYNCCWDQCQACFNSSPDLADHIRSIHVDGQRGGVVFVCLWKGCKVYNTPSTS QSWLQRHMLTHSGDKPFKCVVGGCNAFASQGGGLARHVPTHFSSQNSSKVSSQPKAKEESPKAGMKNRRKLNKRRR RSLPRPHDFFDAQTLDAIRHRAICFNLSAHIESLGKGHSVVFHSTVIAKRKEDSGKIKLLHWMPELDPDVVWVNESERHQL KTKVVHLSKLPKDTALLLDPNIYRTMPQKRLKR

靶点信息

研究背景	Acts as an accessory subunit for the core Polycomb repressive complex 2 (PRC2, which mediates histone H3K27 (H3K27me3 trimethylation on chromatin leading to transcriptional repression of the affected target gene. Plays a role in nucleosome localization of the PRC2 complex.
------	---

基因ID	121536
基因名	AEBP2
Swiss	Q6ZN18
别名	AEBP2

产品验证



Western blot analysis of AEBP2 expressed in Mouse kidney,Rat heart,Jurkat using AEBP2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn