

ETS2 Rabbit pAb

货号: **AYP16379**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	53kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,mouse lung
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-469 of human ETS2 (P 15036).
序列	MNDFGIKNDQVAPVANSYRGT LKRQPAFD TFDGSLFAVFP SLN EEQTLQEVPTGLDSISHDSANCELPLLT PCSKAVMS QALKATFSGFKKEQRRLGIPKNPWLWSEQVCQWLLWATNEFSLVNVNLQRFGMNGQMLCNLGKERFLELAPDFVGD I LWEHLEQMIKENQEKTEDQYEENSHLTSVPHWINSNTLGFGEQAPYGMQTQNYPKGGLLDSMCPASTPSVLSSEQEFQ MFPKSRLSSSVTYCSVSQDFPGSNLNL TNNSGTPKDHDSPENGADSFESSD SLLQSWNSQSLLD VQRVPSFESFED DCSQSLCLNKPTMSFKDYIQRSDPVEQ GKPVIPAAVLAGFTGSGPIQLWQFLELLSDKSCQS FISWTGDGW EFKLADP DEVARRWVGKRKNPKMNYEKL SRGLRYYYDKNIIHKTS GKRYVYRFVCDLQNL LGFTPEELHAILGVQPD TED

靶点信息

研究背景	This gene encodes a transcription factor which regulates genes involved in development and apoptosis. The encoded protein is also a protooncogene and shown to be involved in regulation of telomerase. A pseudogene of this gene is located on the X chromosome. Alternative splicing results in multiple transcript variants.
基因ID	2114
基因名	ETS2
Swiss	P15036
别名	ETS2;ETS2IT1

产品验证



Western blot analysis of ETS2 expressed in A-431, mouse lung using ETS2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn