

SMYD4 Rabbit pAb

货号: **AYP16386**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:100
理论分子量	89kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse brain,Mouse ovary,Rat kidney,Rat brain
细胞定位	cytoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human SMYD4 (NP_443160.2).
序列	MDLPVDEWKS YLLQKWASLPTS VQVTISTAETLRDIFLHSSLLQPEDELFLKRLSKGYLVGKDS DAPLFYREEGNKKFQE KDYTGAAVLYSKGVSHSRPNTEDMSLCHANRSAALFHLGQYETCLKDINRAQTHGYPERLQPKIMLRKA ECLVALGRLQEA SQTISDLERNFTATPALADVLPQTLQRNLHRLKMKMQEKDSLTESFPAALAKTLEDAALREENEQLSNASSSIGLCVDPLKG RCLVATKDILPGELLV

靶点信息

研究背景	Plays a critical role in cardiac development. Acts as a key epigenetic regulator of gene expression during cardiac development via its dual activities as a methyltransferase and negative regulator of HDAC1 (By similarity).
------	--

基因ID	114826
基因名	SMYD4
Swiss	Q8IYR2
别名	SMYD4;ZMYND21

产品验证



Western blot analysis of SMYD4 expressed in Mouse kidney,Rat kidney,Mouse brain using SMYD4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 12 0s.

实验步骤

访问官网浏览详情: www.ablybio.cn