

# TP53AIP1 Rabbit pAb

货号: **AYP16424**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	8kDa/11kDa/12kDa
实测分子量	13kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse liver,Rat liver
细胞定位	Mitochondrion
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-108 of human TP53AIP1 (NP_001238893.1).
序列	MGSSSEASFRSAQASCSGARRQGLGRGDQNLVMPNGRAQTHTPGWVSPCSENRDGLLPATAPGRLCSHRGADIPS FQTHQDPVTASGSSELHADCPQFRALDRAGN

## 靶点信息

研究背景	This gene is specifically expressed in the thymus, and encodes a protein that is localized to the mitochondrion. The expression of this gene is inducible by p53, and it is thought to play an important role in mediating p53-dependent apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.
基因ID	63970

基因名	TP53AIP1
Swiss	Q9HCN2
别名	TP53AIP1;P53AIP1

## 产品验证



Western blot analysis of TP53AIP1 expressed in Mouse liver, Rat liver using TP53AIP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)