

# SOX14 Rabbit pAb

货号: **AYP16426**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	26kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 91-240 of human SOX14 (NP_004180.1).
序列	RYVFPLPYLGDTDPLKAAGLPVGASDGLLSAPEKARAFPPASAPYSLLDPAQFSSSAIQKMGEVPHTLATGALPYASTLGY QNGAFGSLSCPSQHTHTHPSTNPGYVVPNCNCTAWSASTLQPPVAYILFPGMTKTGIDPYSSAHATAM

## 靶点信息

研究背景	This intronless gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. Mutations in this gene are suggested to be responsible for the limb defects associated with blepharophimosis, ptosis, epicanthus inversus syndrome (BPES) and Mobius syndrome.
基因ID	8403

基因名	SOX14
Swiss	O95416
别名	SOX14;SOX28

## 产品验证



Western blot analysis of SOX14 expressed in Mouse brain,Rat liver using SOX14 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)