

— ABLYBIO, Help Your Research



# ELMO3 Rabbit pAb

货号: **AYP16442**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	68kDa/81kDa/87kDa
实测分子量	81kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	22Rv1,SW480,HT-29,Mouse brain,Mouse liver,Mouse kidney,Rat lung
细胞定位	Cytoplasm
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 504-773 of human ELM O3 (NP_078988.2).
------	--

## 靶点信息

研究背景	The protein encoded by this gene is similar to a <i>C. elegans</i> protein that functions in phagocytosis of apoptotic cells and in cell migration. Other members of this small family of engulfment and cell motility (ELMO) proteins have been shown to interact with the dedicator of cyto-kinesis 1 protein to promote phagocytosis and effect cell shape changes.
基因ID	79767
基因名	ELMO3
Swiss	Q96BJ8 ( <a href="https://www.uniprot.org/uniprotkb/Q96BJ8/entry">https://www.uniprot.org/uniprotkb/Q96BJ8/entry</a> )
别名	ELMO3,CED-12,CED12,ELMO-3,ELMO3 Rabbit pAb

## 产品验证



Western blot analysis of ELMO3 expressed in Mouse brain, Rat lung, 22Rv1 using ELMO3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))